


iRhythm™ INTEGRATION SERVICE

iRhythm will offer Electronic Health Record (EHR) integration designed to streamline your ZIO® Service workflows:

- Centralize workflows around your EHR system by enabling ordering and receiving of the ZIO Service through your EHR
- Minimize chances for human error by eliminating double data entry into ZIOreports.com and your EHR
- Speed up ZIO Patch Report publishing and association with the patient medical record in your EHR

iRhythm Integration Service is designed to fit smoothly within your institution's existing workflow in five convenient steps:



Key advantages of implementing the iRhythm service:

Efficient Workflow and Reduced Administrative Burden

- » ZIO Service ordering and automatic patient registration from within your EHR
- » ZIO Patch Reports published directly to your EHR and automatically associated to patient medical records in a timely fashion

Greater Accuracy in Patient Data Management

- » Latest revision of ZIO Patch Report automatically published to your EHR


Ease of Integration

- » Onsite hardware or software not required
- » Fully dedicated iRhythm technical project manager

Compatibility and Standardization

- » Compatible with any EHR system that processes HL7 messages

EHR INTEGRATION DATA FLOW¹




¹Integration Service will also be available with the option for physicians to interpret ZIO Patch Reports in ZIOreports.com to preserve current workflow. Physician interpreted reports will then be posted to your EHR system.

System Requirements

- HL7 Version 2.x or higher
- Secure communication interface (e.g., SFTP, VPN)
- Support for base64 encoded PDF file format

Implementation Timeline²



² Ideal timeline represented; actual schedule to be agreed upon with your organization's IT resources.

D100A4003.01